




CarbonFore

No More Stress on the Nose Bridge
Compatible with Your Glasses

*Unique Forehead Frame Design
for Stress-Free Comfort*





Unique Forehead Frame Design for Stress-Free Comfort


The Xpedent CarbonFore Loupes are revolutionising comfort and functionality for dental professionals. Designed with a unique forehead frame, these loupes transfer the weight away from your nose, eliminating pressure and providing unparalleled comfort during long procedures—all while allowing you to wear your own glasses.



CarbonFore Keplerian Ergo Loupes

Improve Your Posture

Enhance Your Well-Being



Model	Magnification	Working Distance	Pupillary Distance Adjustment	Depth Of Field	Field Of View	Weight
FE3.5	3.5X	280-380mm	50-80mm	100mm	70-90mm	79g
FE4.5	4.5X	360-460mm		100mm	50-80mm	81g
FE5.5	5.5X	440-540mm		90mm	50-75mm	87g
FE6.5	6.5X	500-600mm		80mm	45-70mm	85g

Why Choose CarbonFore Loupes?

No Stress for the Nose

Crafted with an ultra-light alloy structure, the CarbonFore Loupes minimise weight and eliminate nose pressure, ensuring unparalleled comfort during even the longest procedures. Experience the freedom of a loupe that feels almost weightless yet delivers top-tier performance.

Unique Forehead Frame Design

Unlike traditional loupes, the CarbonFore Loupes feature a specially engineered forehead frame that redistributes the weight evenly across your forehead. This innovation minimises strain on your nose, ensuring a stress-free and comfortable experience even during extended use.





Why Choose CarbonFore Loupes?

Works with Your Glasses

Designed for ultimate compatibility, the CarbonFore Loupes clip securely onto your own prescription glasses, allowing you to enjoy both optical precision and the familiarity of your own lenses.

Ergonomic and Flexible

The Flip-Up Function provides easy transitions between magnified and natural vision, ensuring unobstructed patient interaction and a seamless workflow.

Why Choose CarbonFore Loupes?

Magnification Options

Available in four sizes—3.5x, 4.5x, 5.5x, and 6.5x—the CarbonFore Loupes offer the perfect magnification level for every dental procedure.

Advanced Optics

Featuring Schott Glass and Achromatic Optical Path Design, these loupes deliver crisp, high-contrast images with exceptional depth of field and true color fidelity, enabling precise visualisation.

Ultra-Lightweight Construction

Crafted from ultra-light alloy, the CarbonFore Loupes combine a sleek metallic look with outstanding durability and minimal weight for maximum comfort.





Reduces Eye Fatigue from Long Surgeries

Independent Adjustment of Left and
Right Interpupillary Distance
Easy Image Alignment



Experience Stress-Free Dentistry

The Xpedent CarbonFore Loupes redefine comfort and performance with their innovative forehead frame design and ergonomic features. Say goodbye to nose strain and hello to a stress-free, efficient, and precise dental experience.



CarbonFore Keplerian Loupes

Model	Magnification	Working Distance	Depth Of Field	Field Of View	Weight
F4.0	4.0X	280-380mm	200mm	150-170mm	84g
F5.0	5.0X	360-460mm		130-150mm	86g
F6.0	6.0X	500-600mm		110-130mm	94g



CarbonFore Galilean Loupes

Model	Magnification	Working Distance	Pupillary Distance Adjustment	Depth Of Field	Field Of View	Weight
F2.5	2.5X	300-580mm	50-80mm	200mm	150-170mm	71g
F3.0	3.0X			200mm	130-150mm	77g
F3.5	3.5X			200mm	110-130mm	81g

Xpedent CarbonFore Loupes

No Pressure Just Precision



xpedent_uk



Xpedent International



Xpedent



jenny@xpedent.com

www.xpedent.com

info@xpedent.co.uk