



Fire fly

Cordless Headlight



55000 LUX

Fire fly provides stronger light intensity to enable work with higher magnifications.



100mm



CRI>90%

Lightweight

Optimised weight distribution, ergonomically designed to rest lightly on your loupes





Cordless

Freedom of movement increases comfort and efficiency

Compatible

Can be mounted on most loupes models



User-friendly controls

Simple operation with tap controls
for ON/OFF and light intensity



600mAh

High-Capacity Lithium Battery for Extended Battery Life



touching area

- Touch once for the brightest setting.
- Touch twice for medium brightness.
- Touch three times for the lowest brightness.
- Touch four times to turn off the power.

Swift battery replacement

Enabled by magnetic mount.






Angle Adjustment

The headlamp can be adjusted up and down



Parameter

Intensity	55K LUX
Weight	27 g
Run time	2 hrs
Charge time	1.25hrs
Batteries	3 pcs
Colour temp'	5,750K
Yellow filter	External
Colour	 BLUE  GOLD  GRAY





Follow up



xpedent_uk



Xpedent International



Xpedent



jenny@xpedent.com

www.xpedent.com **info@xpedent.com**



Xpedent Head Office

Xpedent UK Ltd

88 Boundary Road,
Hove, East Sussex,
BN3 7GA
England 06128867



Sales Office

Xpedent Austria GmbH

Millennium Park 9
6890 Lustenau
Austria



Manufacturer

Guilin Yikeshi Medical
Instrument Co., Ltd

D08 Building, High Tech
Industrial Park, Chaoyang Road
Qixing District Guilin, 541004
China